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Hersteller:
Dr. Jean Bausch GmbH & Co. KG
Oskar-Schindler-Str. 4
50769 Köln - Germany
Telefon: +49 221 70936-0
Fax: +49 221 70936-66
E-Mail: info@bauschedental.de
Internet: www.bauschedental.de

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Instructions for use as at: 23-05-14

1 General description and intended purpose of the medical product

Sterilisable, flexible blue silicone sleeves, 45 mm long, with an external diameter of Ø 4.5 mm.

Indication

Arti-Grip sleeves are used to improve the clamping effect of articulating-paper forceps (Miller forceps). The silicone sleeves, pre-cut in 45 mm lengths, are fitted over one or both clamping surfaces of Miller forceps. The slightly adhesive surface of the silicone sleeve enhances the clamping effect to ensure that thin occlusion foils and shimstock foils are gripped reliably. Arti-Grip sleeves can be used by dentists in the mouth of the patient, on models and as well by dental technicians in the laboratory. In situations where special hygiene requirements are indicated for a patient, Arti-Grip sleeves can be sterilized prior to use.

Contraindication

No known contraindications

Side-effects

Possible side-effects may be allergic reactions. Experiences with the products and clinical evaluation indicate that no negative reports could be verified in relation to an allergic reaction. There has not been a single justified complaint in relation to an allergic reaction. If these products are used in conjunction with new, as yet unknown products or materials, side-effects or interactions may occur.

2 Information

- The instrument can be sterilised using steam (moist heat, 134°C).
- The Arti-Grip sleeves must always be removed from the clamping surfaces of the forceps and then cleaned, disinfected and sterilised before use.
- Disposal: Before disposal the Arti-Grip sleeves must be reconditioned (sterilised) or disposed of with the usual contaminated waste from the dental practice.

3 Pack contents

- 20 silicone sleeves (non-sterilised)
- Instructions for use

4 Preparation

Arti-Grip sleeves must be cleaned, disinfected and sterilised before initial use and before each further use. Please also note the detailed guidelines regarding reconditioning at Point 6 "Reconditioning".

5 Application

- Arti-Grip sleeves must be removed from the sterile packaging following the usual practice hygiene measures (use of disposable gloves).
- The silicone sleeves should be wetted with medical alcohol (ethanol or isopropyl alcohol min. 70%) beforehand. Once wetted, the sleeves are easier to slip over the clamping surfaces of the forceps.
- A strip of articulating material (paper or foil) is then clamped between the clamping surfaces of the forceps.
- The articulating material is retained between the relevant teeth of the upper and lower jaw with the holder (buccal).
- Due to the increased clamping force a small strip of shimstock foil (BK 38) can also be fixed in position between the two tips of the clamping surfaces
- Static and dynamic occlusion can be tested.
- The articulating material should be removed from the mouth and disposed of with the usual contaminated waste from the dental practice.
- The holder and silicone sleeves are then prepared for reconditioning.

6 Reprocessing

The following should always be noted: "Hygiene requirements during the preparation of medical products recommendation of the Commission for Hospital Hygiene and Infection Prevention (KRINKO) at the Robert Koch Institute (RKI) and the Federal Institute for Drugs and Medical Devices (BfArM)", Bundesgesundheitsblatt (Federal Health Gazette) 2012 • 55:1244–1310.

The instruments and Arti-Grip silicone sleeves must be reconditioned after each use:

- Remove the silicone sleeves from the clamping surfaces. Thoroughly hand-clean the parts with correspondingly small brushes with firm bristles. The inside of the Arti-Grip silicone sleeves must also be cleaned using suitable instruments.
- Information: The parts should be cleaned in a water bath below the water surface without any other cleaning additives to achieve both adequate cleaning of the instruments yet prevent protein fixation and protect the surrounding area against contamination with water spray.
- Rinse the instrument parts intermittently using water (drinking water quality minimum).
- Place the parts in a suitable cleaning and disinfectant solution routinely used in the practice. Examples:
 - o Becht Bechtol Futura
 - o Dürr Dental ID 213 instrument disinfection
 - o Pluradent Pluline instrument solution
 - o Schülke & Mayr Lysetol V
- o See also "List of disinfectants and disinfection procedures certified and recognised by the Robert-Koch Institute" or the VAH list of disinfectants.
- o Note: Cleaning and disinfectant agents containing formaldehyde may only be used after adequate cleaning to prevent protein fixation.
- o Note: Strictly follow the instructions provided by the manufacturer of the cleaning agent or disinfectant. It is particularly important to observe the concentrations to be used and reaction times!
- o Information: For mechanical cleaning, the manufacturer's instructions for use must be strictly observed! With haemostatic forceps, ensure that the clamping surfaces are exposed and not covered!
- o Finally, rinse the instrument parts with water (drinking water quality minimum, recommended: demineralised water with microbiological quality corresponding to drinking water).
- Dry.
- Visual check for contamination. Should residual contamination persist, the entire cleaning procedure with all steps must be repeated (cleaning, intermittent rinsing, disinfection, final rinsing and drying).
- The instrument and all parts must be free of any residue and dry before further preparation.
- Pack and seal all parts of the instrument in suitable, fitting clear sterile packaging (sterilisation bags) in compliance with EN 868-5. Observe the instructions of the manufacturer of the sterilisation bag and sealing machine and the current standard requirements.

- Sterilisation must be completed in a validated procedure using moist heat in an autoclave. Observe the instructions for use of the autoclave manufacturer.
- Sterilise the instrument using moist heat (saturated steam) using a pre-vacuum procedure for **5 minutes at 134°C**.
- After completion of the sterilisation, the sterile packaging must be marked with the date up to when the instrument can be regarded as sterile assuming the packaging remains undamaged.
- The recommended storage period for sterile medical products is described in Standard DIN 58953-8 and depends on external influences and effects during storage, transport and handling.
- Keep the packaged, sterile instrument protected in a dry, clean place.

7 Symbols

Sterilise with moist heat (autoclave) at up to 135°C

Bausch Articulating Papers, Inc.
12 Murphy Drive, Unit 4
Nashua, NH 03062
United States of America
Phone: (603) 883-2155
Phone: 1-888-6-BAUSCH
Fax: (603) 883-0606
Internet: www.bauschedental.com
E-Mail: info@bauschedental.com

Bausch Articulating Papers (Australia) Pty. Ltd
G.P.O. Box 3733
Sydney NSW 2001, Australia
Tel: +61-2-9345-1945
Fax: +61-2-9345-1955
Internet: www.bauschedental.com.au
E-Mail: info@bauschedental.com.au

Bausch Articulating Papers Japan K. K.
2nd Floor, 1-4-2 Jonan,
Ikedashi, Osaka Japan 563-0025
T: +81 72-737-9501
F: +81 72-737-9502
E-Mail: info@bauschedental.jp
Web: www.bauschedental.jp

Bausch Importação de Materiais Odontológicos Ltda. EPP
Rua Paulo Eduardo Xavier de Toledo, 379 salas 8 e 9 Lt. 3 Q. 30.
13304-240 Itu-SP, Brasil
Tel: +55-11-4813 2806
Fax: +55-11-4813 2806
E-Mail: vendas@bausch.net
Internet: www.bausch.net

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